DSR 130

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Dear users,

we congratulate on the buying of a technically matured and the best available technology product.

This instructions helps you and the careers as a basis for sure, safe use of the Shower chair and a sensible use of the product possibilities.

The instructions describes the preparation, service, use, servicing and cleaning of the product.

This instructions is to be perused carefully and must be always accessible to all users of the product.

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Standing insecurity / standing inability with use decrease / of missing functional ability of the lower extremities or other illnesses who are affecting the stability (e.g., strength loss with evil illnesses, chronic, advanced heart circulatory illnesses and lung illnesses) infers with a body position weakness, so that an independent use of the shower is still allowed.

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This product is suitable for reapplication. With the processing the standards of the manufacturer must be considered. This standards are made available on inquiry.

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- Read this instructions carefully.
- Before use all components are to be checked for damages and firm seat.
- While seating, rising and while folding up or engaging the seat plate it's possible to pinch clothing or body parts
- The steadiness of the shower chair is clearly lower as with a normal chairs, therefore avoid lateral lean out – Tilting danger !
- The Shower chair DSF 130 serves excluding the use by a person to max. 130 kg on the seat.

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
The shower chair is preassembled, the backrest must be only be latched in and the 4 foots must be assembled (see next chapter).

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To adjust the footrests in the height press the brass pin. Then the footrests can be shifted in 5 steps in the height.

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- 1 Product name
- 2 Article number
- 3 Serial number and date of production (month/Year)
- 4 Attention: See operating instructions

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